

ChIP assay - Swaroop Lab, 2013

ChIP protocol (Sooyoung) - For retina tissues (Check: one retina amount of gDNA:)

Generally, use DNA, 5, 10, or 20 ug, but some peoples are using 200 ug also.

H3, H4, 2ug, 2hr, 5 ul PCR - result good

Generally, other antibody - DNA 20 ug, 2 ug antibody, would be OK, but need optimization

Prior to Day1- option

Coating protein A beads

Take enough dynabeads slurry (30ul-50ul each sample) (Invitrogen #100-01D)

Separate beads and buffer with Magnetic Separation Stand (promega # Z5332)

Remove the sup. with loading tip

Wash three times with 1 mL IP dilution buffer

Resuspend beads in 1mL pre-blocking buffer

Incubate at 4°C, greater than 2hr or O/N

Separate beads and buffer with Magnetic Separation Stand and remove the sup.

Wash twice with 1 mL IP dilution buffer

Add IP dilution buffer back to original volume

(Comments: This protocol is modified specifically for ChIP-Seq. Many protocols suggested that dynabeads don't need blocking. For separation of beads and sup, I first briefly centrifuge it and then put the tube to the magnetic stand. Remove the sup. by loading tips)

DAY1

Cell lysis and chromatin DNA extraction

- 1) Dissect retinas from eye balls, and wash with Cold PBS (X1)
- 2) Put retinas into 900 ul RT PBS
- 3) Add 100 ul formaldehyde (from 10% stock) and mix immediately, Incubate and rotate at RT for 15 min
For histone antibody, 10 min would be proper, but for other antibody the time could be increased upto 30 min. (Check, 10 min, 15 min and 30 min)
- 4) Add 143 ul **glycine** (from 1M Stock, final 125 mM) and mix immediately and spin down, remove sup
- 5) Add **0.125 M glycine** to soak retina, rotate 5 min at RT

Keep 4 °C (on ice) from here

- 6) Rinse retinas with 1ml Cold 1xPBS twice, centri at 4°C
- 7) Add **600 ul Cell lysis buffer** and keep on ice for 10 -20 min and Pipet up and down every 5-10 min until no clear tissue can be visible (or use homogenizer)
- 8) Centri 2500g (3500 rpm), 5 min in Cold centrifuge and remove the sup.
- 9) Resus pellet in at least 10 vol (600 ul) with cell lysis buffer and Rinse once and centri 2000 g (3500 rpm) again, remove the sup.
- 10) Re-suspend the nuclear pellet (white) in 200 ul **nuclei lysis buffer** and pipet up and down to disrupt clumps
- 11) On ice 10-20 min (lay the tube on ice to avoid precipitation of SDS in the nuclei lysis buffer)
- 12) Add 400 ul **IP dilution buffer**+protein inhibitors

Possible to freeze at -80°C

Sonication

- 13) Sonicate –Covaris setting (*check: and gel loading*)
Sonication setting (___ pulses of ___ sec, X___, interval ____,)
Take 10 ul for agarose gel electrophoresis to check the size of fragmented DAN after decross-linking
Average fragment size should be between 100 bp-600 bp.
- 14) Centri at 14000 rpm 15 min (5 min, or 10 min, OK)
- 15) Take the supernatant and centri again at 14000 rpm for 15 min (5 min or 10 min, OK even once centri would be OK). Save 10ul for input DNA (1/10 of working)

Possible to freeze at -80°C

Antibody binding

- 16) Preclear the lysate by incubating by constant rotation with dynabeads (Invitrogen, 25 ul) for >2 hours in the cold room.

ChIP assay - Swaroop Lab, 2013

- 17) Separate beads and buffer with Magnetic separation stand.
- 18) Take the sup.
- 19) Add antibody 5 ug to the sample (Keep the no-antibody control sample at 4°C without adding any antibody), rotate O/N at 4°C.

DAY 2

Dynabeads binding, washing, elution, and decross-linking

- 20) Add 30 ul blocked dynabeads to each sample (also to no-antibody control tube)
- 21) Incubate on a rotating wheel at RT 60-90min (or cold room, 6 hrs)
- 22) Separate beads and buffer with Magnetic separation stand.
- 23) Remove sup (Beads from no-antibody tube could serve as "no Ab" negative control after washing and elution)
- 24) Wash the beads twice with 1 ml 1X dialysis buffer
Add 1 ml buffer, rotate 15 min at RT, separate beads and buffer with Magnetic separation stand, remove the sup. As much as possible with loading tip
- 25) Wash the beads twice with 1ml IP wash buffer (same as above)
(comment: Washing step could be held at 4°C with longer time, if non-specific binding still occur, wash 3 times for each buffer followed by 1-3 times TE wash)

Elution and decross-linking

- 26) Add 150 ul elution buffer to the beads
- 27) Incubate in 65°C water bath for 10-15 min (vortex every 2-3 min)
- 28) Separate beads and buffer with Magnetic separation stand
- 29) Transfer the sup. To a new tube.
- 30) Repeat elution steps again and combine both elution (300ul)
- 31) Add 30 ul 3M NaCl (final 250mM) and 1 u 10mg/ml RNaseA (for input sample, add right amount of NaCl and RNase A)
- 32) Incubate 65°C, 4-5 hr or O/N (input control samples also need to be decross-linking at the time)
- 33) Purify ChIP DNA and input DNA with Qiagen PCR purification kit
- 34) Do qPCR with ChIP sample, no-Ab control and total input

Recipes for solution and material information

1 M glycine (FW = 75.07)

3.75 g / 50 mL dH₂O

filtrate and store at RT

1% formaldehyde in 1x PBS (freshly prepared)

Total 25ml

37% formaldehyde 676 uL

1x PBS 24.324ul

0.125 M glycine in 1x PBS (freshly prepared)

10 mL 30 mL

1 M glycine 1.25 mL 3.75 mL

1x PBS 8.75 mL 26.25 mL

Proteinase inhibitors

(final conc.)

(per 1 mL solution/buffer)

100 uM PMSF 1 uL of 100 mM

100 ug/mL benzamideine 10 uL of 10 mg/mL

5 ug/mL leupeptin/pepstatin 1 uL of 5 mg/mL

Cell lysis buffer (at 4°C) (40 mL)

10 mM Tris-Cl (pH 8,1) 400 uL of 1 M (Using Tris Base)

10 mM NaCl 133 uL of 3 M

0.5% NP-40 2 mL of 10%

proteinase inhibitors (freshly add)

dH₂O 37.47 mL

ChIP assay - Swaroop Lab, 2013

Nuclei lysis buffer (20 mL)
50 mM Tris-Cl (pH 8.1) 1 mL of 1 M
10 mM EDTA 400 uL of 0.5 M
1% SDS (Sodium Dodecyl Sulfate) 2 mL of 10%
proteinase inhibitors (freshly add)
dH₂O 16.6 mL

IP dilution buffer (at 4°C) (40 mL)
16.7 mM Tris-Cl (pH 8.1) 668 uL of 1 M
167 mM NaCl 2227 uL of 3 M
1.2 mM EDTA 96 uL of 0.5 M
1.1% Triton X-100 2.2 mL of 20%
0.01% SDS 40 uL of 10%
proteinase inhibitors (freshly add)
dH₂O 34.77 mL

Pre-Blocking buffer for dynabeads: 1mL
glycogen 10uL (20mg/mL),
BSA 40uL (5mg/mL)
Yeast-tRNA 10uL (20mg/mL)
in 1mL IP dilution buffer

500 mM NaHCO₃ (filtrate, store at RT)

NaHCO₃ (FW = 84.01) 4.2 g
H₂O up to 100 mL

1x dialysis buffer (store at 4°C) (100 mL)
50 mM Tris-Cl (pH 8.1) 5 mL of 1 M
2 mM EDTA 400 uL of 0.5 M

ChIP assay - Swaroop Lab, 2013

0.2% sarkosyl (omit for monoclonal Ab) 1 mL of 20%

The other name is N-lauroylsarcosine solution (SIGMA L7414-10ML)

dH₂O 92.6 mL

IP wash buffer (stored at 4°C) (100 mL)

100 mM Tris-Cl (pH 9.0) (pH 8.0 for monoclonal Ab) 10 mL of 1 M

500 mM LiCl 10 mL of 5M

1% NP-40 10 mL of 10%

1% deoxycholic acid (Sodium deoxycholate Sigma D6750-25G) 5 mL of 20%

dH₂O 65 mL

Elution buffer (40 mL)

50 mM NaHCO₃ 4 mL of 500 mM

1% SDS 4 mL of 10%

dH₂O 32 mL of dH₂O

5x PK buffer

50 mM Tris-Cl (pH 7.5) 5 mL of 1 M

25 mM EDTA 5 mL of 0.5 M

1.25% SDS 12.5 mL of 10%

dH₂O 77.5 mL

Glass beads, acid-washed, Sigma G1277-100G

10 mg/mL proteinase K

10 mg/mL glycogen

3 M NaCl

Phenol (saturated with TE)

ChIP assay - Swaroop Lab, 2013

CHCl_3 (+1/25 vol. Iso-amyl alcohol)

EtOH (-20°C)